

MADE IN FRANCE 

1 – Device identification :

CE Marking	Index of classification	Notified body	PHTALATES LATEX	
CE 0459	Ila	G-MED		

The pressure cuffs are medical devices used to pressurize, under control of a manometer, a pocket of product to be infused or transfused quickly.

The specific design of the pressure cuff can be wrapped around the liquid bag and thereby makes it very easy implementation and cleaning. The pressure is controlled by manometer. Depression is controlled by a valve.

2 – References :

Reference	Designation	Size
M30500	Transparent polyurethane cuff	500 cc
M31000	Transparent polyurethane cuff	1000 cc
M34000	Transparent polyurethane cuff	3000 cc



3 – Technicals characteristics :

- Pressure cuff : TPU
- Manometer : ABS + glass plexiglass
- Manometer protector : Silicone
- Valve bulb : Chrome plated brass
- Bulb : PVC phtalates free

4 – Operating instructions :

- Open the pressure cuff.
- Place the liquid bag to the center of the pressure cuff, output tubes oriented to bulb.
- Insert the fastening system of pressure cuff in the opening of the liquid bag before insert it into the opening made on the pressure cuff.
- Close the pressure cuff by squeezing the movable portion around the liquid bag.
- Using the bulb, breathe the pressure required to maintain the cuff around the liquid bag, suspend all on a drip stand.
- As to the flow of product irrigated, the pressure displayed by the manometer falls gradually : it is necessary to restore order to ensure complete emptying of the liquid.

5 - Maintenance :

- Calibration control of the manometer : every 12 months
- Visual inspection : before use the medical device.

6 - Cleaning :

- Clean the surface of the complete device.
- Set can be decontaminated by spraying disinfectant solution

The cleaning is done with :

- SURFA'SAFE Premium
- Or
- ANIOXY-SPRAY WS

Disinfection is done with :

- ANIOSPRAY SURF 29

7 - Storage conditions :

Type of packaging	Storage	Temperature	Humidity	Atmospheric pressure
Original packaging	Ventilated area	-10° à 40 ° c	30 to 40 %	500 to 1060 hPa

8 – Warranty :

This warranty provides assurance for the customer who purchases a D & D product that should the product fail to function to D & D published specifications during the term of this warranty, will either replace or repair.

The guarantee period is 2 (two) years from the date of purchase.

The product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.

9 - Warning :

Do not use the device near a heat source.

This use could cause device failure and a risk of complications for the patient.