

尺寸：148x210mm  
材质：128g双铜  
工艺：双面印刷+四色

feellife

The third way of dosing



Hidden medication cup

Eliminate the fear of nebulization for children

Cute cartoon + Breathing light

Attract children's attention

Nearly silent operation, multi-angle nebulization

Suitable for using when children are sleeping.

Air Angel

Make children enjoy nebulization

Small and Uniform particles

Can be absorbed more easily for children

Musical chip built inside

Compatible for downloading children preferred music from computer, making it a pleasant time when use.

Technical parameters

Model: Air Angel

Power Supply: Internal power supply

Power Consumption: <3.0w

Nebulization Rate: >0.3ml/min

MMAD:  $2.5 \pm 30\% \mu\text{m}$

Medication Cup Capacity: 6ml

Device Dimensions: 53(L) \* 52(W) \* 140(H)mm

Device Weight: 107g



Deposition site: reach to bronchia and alveoli

Authoritative test report of proportion of particle size

The medium position particle size of the nebulizer is in the temperature of  $25 \pm 3^\circ\text{C}$  and the humidity of  $60 \pm 5\% \text{RH}$  with the 0.9% physiological saline. Under these conditions the particle size distribution curve is as follows:



Notes: The horizontal axis is the particle size which shows Logarithmic distribution. The vertical axis on the left is the cumulative percentage of the volume, corresponding to the upward trend of the graph. The vertical axis on the right is the volume percentage of a certain interval, corresponding to the histogram or undulating curve.

\* Application: To make the liquid medication nebulized so that the patients can inhale directly.



FEELLIFE HEALTH INC.

Factory: Building A, NO. 9 Furong Road, Songgang Street, Bao'an District, Shenzhen 518104

Website: [www.feellife.com](http://www.feellife.com)